




Scan to read an in-depth version of this poster.

Einar Halldórsson¹, Emil G. P. V. Stender, Ph.D.¹

(1) Fida Biosystems, Copenhagen (Søborg), Denmark, www.fidabio.com, e-mail: einar@fidabio.com; egps@fidabio.com

Method introduction

Flow Induced Dispersion Analysis (FIDA) is a first-principle technique used to measure the hydrodynamic radius of molecules. FIDA is highly sample-efficient, requiring only 40 nL of sample per measurement, and operates without any buffer constraints. Recently, specialized methods to study amyloid fibrils and protein condensate have been developed with the Fida instruments (MTIS¹, Capflex², TDIPS³).




Studying Fibril Kinetics & Interactions quantitatively with FIDA

1. FIDA enables you to **detect and measure** the size of monomers, small oligomers and fibrils with out any size bias. From the hydrodynamic radius determined by FIDA, the **length of the fibrils** can be calculated.
2. It monitors the development (size increase) of fibrils over time, allowing you to **capture fibril growth kinetics**.
3. It allows you to monitor fibril **growth inhibition**
4. It enables you to work directly on patient-derived fibrils e.g. cerebral spinal fluid, with only **nanolitres to microlitres of sample required**.

Fibril polymorphism with MTIS

1. MTIS performs an **automated ΔG** comparison for diverse polymorphs.
2. It seamlessly integrates the precision of separation-based techniques with the **speed and efficiency** of bulk methods.
3. MTIS quantifies **both monomer concentration and total concentration in the same assay**.



MTIS compared to other methods to study thermodynamic stability




Figure 2: a-b) Thermodynamic stability data from MTIS of α Syn fibrils formed under different conditions compared to AFM1 and ΔG values calculated and compared⁶. c) Thermodynamic stability of three different amyloid fibrils measured with different techniques and compared to MTIS (FIDA)¹.

Studying protein condensates with Capflex and TDIPS in FIDA

1. Capflex allows you to measure dilute phase concentration, along with the **quantification of number and size** of the condensates. It can also account for the aging of condensates.
2. TDIPS allows you to screen conditions for condensate formation using **nL of your sample**. It also allows you to record phase diagrams label free with nL of your sample.




Figure 3: a) Phase diagrams of α Syn variants recorded with TDIPS, confirmed by confocal microscopy. b-d) Dilute phase concentration of α Syn variants measured with Capflex and confirmed by confocal microscopy⁷.




Figure 1: a) Affinity determination of antibodies against fibril oligomers⁴. b) Binding to full length fibrils⁵. c) Monitoring fibril growth and fibril growth inhibition by FIDA using Fibril Paint⁵.

References

1. Farzadfard, A. et al. Thermodynamic characterization of amyloid polymorphism by microfluidic transient incomplete separation. *Chem Sci* 15, 2528–2544 (2024).
2. Stender, E.G.P., Ray, S., Norrild, R.K. et al. Capillary flow experiments for thermodynamic and kinetic characterization of protein liquid-liquid phase separation. *Nat Commun* 12, 7289 (2021). <https://doi.org/10.1038/s41467-021-27433-y>
3. R. K. Norrild, T. O. Mason, L. Boyens-Thiele, S. Ray, J. B. Mortensen, A. W. Fritsch, J. M. Iglesias-Artola, L. K. Klausen, E. G. P. Stender, H. Jensen, A. K. Buell. *Angew. Chem. Int. Ed.* 2024, 63, e202404018. <https://doi.org/10.1002/anie.202404018>
4. Nielsen, J. et al. Molecular properties and diagnostic potential of monoclonal antibodies targeting cytotoxic alpha-synuclein oligomers. *NPJ Parkinsons Dis* 10, (2024).
5. Aragonès Pedrola, J. et al. Fibril Paint To Detect Amyloids and Determine Fibril Length. [doi:10.1101/2023.10.13.562220](https://doi.org/10.1101/2023.10.13.562220).
6. Fricke, C. et al. Thermodynamic Stability Modulates Chaperone-Mediated Disaggregation of Alpha-Synuclein Fibrils. *BioRxiv* 2024.12.19.629136; [doi: https://doi.org/10.1101/2024.12.19.629136](https://doi.org/10.1101/2024.12.19.629136).
7. Ray, S. et al. Divergent effects of pathological Alpha-synuclein truncations and mutations on phase separation. *BioRxiv* 2024.11.18.624073; [doi: https://doi.org/10.1101/2024.11.18.624073](https://doi.org/10.1101/2024.11.18.624073)
8. Princen, K. et al. Pharmacological modulation of septins restores calcium homeostasis and is neuroprotective in models of Alzheimer's disease. *Science*. 2024;384:eadd6260. [doi:10.1126/science.add6260](https://doi.org/10.1126/science.add6260).



Scan for more publications based on FIDA measurements