

### MEASURE



AAV Size



Titer



Polydispersity Index, PDI



Aggregation

# 40nL

Ultra-low  
Sample volume

# 5%

High Sensitivity  
5% size difference  
detectable




Full Automation

2x96 well plates &  
2x50 vials




Exceptional  
Matrix Tolerance

2x96 well plates &  
2x50 vials



### POLYDISPERSITY INDEX, PDI

The Fidabio PDI tool calculates the Polydispersity Index.



### AGGREGATION

The Fidabio Spike counter counts aggregates. Each spike in the fluorescence signal represents one aggregate or one microdrop.